

PROTOCOLS FOR CULTURING *PROCHLOROCOCCUS*

Described here are the protocols we use to grow and cryopreserve *Prochlorococcus*. These protocols work with all strains of *Prochlorococcus* tested to date. We routinely grow *Prochlorococcus* in natural seawater-based (undefined) and artificial seawater-based (defined) media as well as an Instant Ocean-based medium. The natural seawater medium, PRO99, has been described in Moore et al. (2002)¹. The artificial medium for *Prochlorococcus* is a modification of the PCR-Tu₂ recipe from Rippka et al. (2000)² and AMP-1 from Moore et al. (2007)³, and is called AMP-J (Morris et al, 2011)³. We modified the concentrations of metals, EDTA, and bicarbonate to improve the growth of some members of our *Prochlorococcus* stock collection. Trace metals can be toxic to *Prochlorococcus*, so we designed our artificial medium to reduce the chances of metal contamination from impurities in the major salts: we only use ultrapure chemicals, and we reduced the overall metals added to the medium compared to the PRO99 medium. HEPES was replaced with TAPS because HEPES can generate hydrogen peroxide in the light (Morris and Zinser 2013). Measures we take to remove metals from our cultureware (for use with both media types) are described below. Our experience with growing *Prochlorococcus* in large volumes (1L or more) has lead us to modify the culturing procedure, and these modifications are also described. Finally, we describe our methods to cryopreserve *Prochlorococcus* for long term storage and shipment.

CONTENTS:

<u>PROTOCOL</u>	<u>PAGE</u>
PRO99 natural seawater medium	2
AMP-J artificial seawater medium	4
Instant Ocean seawater medium	
Cultureware for <i>Prochlorococcus</i>	7
AMP-J plates	
Growing <i>Prochlorococcus</i> in large volumes (1-10L)	8
Purity test broths	
Cryopreservation methods	10

¹ Moore et al. (2002) Limnol. Oceanogr. 47:989-996.

² Rippka et al. (2000) Int. J. Syst. Evol. Microbio. 50:1833-47

³ Moore et al. (2007) Limnol. Oceanogr. Methods 5:353-362

⁴ Morris et al. (2011) PLoS ONE 6: e16805

⁵ Morris and Zinser (2013) J. Phycol. 49: 1223-1228

I. NATURAL SEAWATER-BASED PRO99 MEDIUM

NOTE: We observe maximal growth with open ocean seawater (i.e. Sargasso seawater), though we can routinely grow *Prochlorococcus* in some coastal seawater as well.

Table 1. Nutrient additions to filtered, autoclaved seawater

<u>Nutrient</u>	<u>Manufacturer/Grade</u>	<u>Primary Stock (M)</u>	<u>Dilution Factor</u>	<u>Final Conc. (μM)</u>
NaH ₂ PO ₄ ·H ₂ O	(1) / ACS	0.025	1:500	50
NH ₄ Cl	(1) / ACS	0.50	1:625	800
Na ₂ EDTA·2 H ₂ O	(2) / 99%	0.012	1:10 ⁴	1.17
FeCl ₃ ·6 H ₂ O	(1) / Analytic	0.012	1:10 ⁴	1.18
ZnSO ₄ ·7 H ₂ O	(3) / >99.5%	0.080	1:10 ⁷	0.008
CoCl ₂ ·6 H ₂ O	(1) / Analytic	0.050	1:10 ⁷	0.005
MnCl ₂ ·4 H ₂ O	(1) / Analytic	0.900	1:10 ⁷	0.090
Na ₂ MoO ₄ ·2 H ₂ O	(2) / ACS	0.030	1:10 ⁷	0.003
Na ₂ SeO ₃	(2) / ~98%	0.100	1:10 ⁷	0.010
NiCl ₂ ·6 H ₂ O	(1) / Analytic	0.100	1:10 ⁷	0.010

Manufacturer Index: (1) Mallinckrodt, (2) Sigma, (3) Fluka

PREPARATION:

1. Prepare a glass filter funnel and flask by cleaning with acid and Milli-Q water, as described on p. 7 (clean the system after every 10 filtrations). Filter raw seawater through a 47mm 0.2 μ m polycarbonate filter (make sure there are no bubbles/creases).

Autoclave seawater in a Teflon bottle (60 minutes for 2L) and allow to cool overnight

2. Prepare primary stocks of NH₄Cl, NaH₂PO₄·H₂O, and trace metals as follows:

A. 0.5M NH₄Cl

- Weigh out 1.34g NH₄Cl using dust-free weigh paper
- Transfer into 50 mL Teflon flask filled with about 40 mL Milli-Q water
- Dissolve NH₄Cl by stirring
- Adjust pH to 7.5
- Transfer to volumetric flask and adjust volume to 50 mL with Milli-Q water
- Filter sterilize with **pre-washed** Millipore Steriflip 0.22 μ M Durapore PVDF membrane system (cat: SE179M6); pre-wash with 50 mL milli-Q water and in laminar flow hood, remove from collection tube; reset tube into filtration system
- In a laminar flow hood, apply sterile tube cap.
- Store sterile stock at 4°C

B. 0.025M NaH₂PO₄·H₂O

- Weigh out 0.15 g NaH₂PO₄·H₂O using dust-free weigh paper
- Transfer into 50 mL Falcon tube filled with about 40 mL Milli-Q water
- Dissolve NaH₂PO₄·H₂O by stirring, inverting tube
- Adjust volume to 50 mL with Milli-Q water
- Filter sterilize with **pre-washed** Millipore Steriflip 0.22 μM Durapore PVDF membrane system (cat: SE179M6); pre-wash with 50 mL milli-Q water and in laminar flow hood, remove from collection tube; reset tube into filtration system
- In a laminar flow hood, apply sterile tube cap.
- Store sterile stock at 4°C

C. Trace Metal Stocks**i. Primary trace metal stocks**

- Using dust-free weigh paper, weigh out:

2.30g	ZnSO ₄ ·7H ₂ O
1.19g	CoCl ₂ ·6H ₂ O
17.81g	MnCl ₂ ·4H ₂ O
0.726g	Na ₂ MoO ₄ ·2H ₂ O
1.73g	Na ₂ SeO ₃
2.38g	NiCl ₂ ·6H ₂ O
- Transfer each trace metal into separate 100mL volumetric flasks containing about 60 mL Milli-Q water
- Dissolve contents by inverting flask several times
- Adjust volume to 100mL mark with Milli-Q water
- Store each stock in a cleaned Teflon bottle at 4°C

ii. Trace metal working stock

- Weigh out 0.435g Na₂EDTA·2 H₂O using dust-free weigh paper
- Transfer into 100mL volumetric flask filled with 60mL Milli-Q water
- Dissolve Na₂EDTA by inverting flask several times
- Weigh out 0.32g FeCl₃·6 H₂O using dust-free weigh paper
- Dissolve FeCl₃ into same volumetric flask by inverting several times
- Individually add and dissolve 100μl of the Primary trace metal stocks
- Adjust volume to 100mL mark with Milli-Q water
- Using a polycarbonate syringe, filter through 0.2μm syringe filter into sterile container in a laminar flow hood
- Store sterile stock at 4°C

3. To make up the media, add following volumes of sterile nutrients and trace metal stock to one liter of the autoclaved seawater. It is important to *dissolve each nutrient sequentially*.

Filtered Seawater	1000.0mL
0.5M NH ₄ Cl	1.6mL
0.025M NaH ₂ PO ₄ ·H ₂ O	2.0mL
Trace metal working stock	100μL

4. Store at room temperature for up to one month.

II. ARTIFICIAL SEAWATER BASED AMP-J MEDIUM

Table 2. Sea Salt and nutrient additions

<u>Nutrient</u>	<u>Manufacturer/Grade</u>	<u>Primary Stock (M)</u>	<u>Dilution Factor</u>	<u>Final Conc.</u>
NaCl	(2) / SigmaUltra	n/a	n/a	481 mM
CaCl ₂ ·2 H ₂ O	(2) / SigmaUltra	n/a	n/a	10 mM
KCl	(2) / SigmaUltra	n/a	n/a	9 mM
Mg SO ₄ ·7 H ₂ O	(2) / SigmaUltra	n/a	n/a	28 mM
MgCl ₂ ·6 H ₂ O	(2) / SigmaUltra	n/a	n/a	27 mM
NaH ₂ PO ₄	(2) / SigmaUltra	0.025	1:500	50 μM
(NH ₄) ₂ SO ₄	(2) / SigmaUltra	0.4	1:1000	400 μM
NaHCO ₃	(2) / SigmaUltra	0.6	1:100	6 mM
TAPS	(2) / SigmaUltra	0.5	1:500	1 mM
Na ₂ EDTA·2H ₂ O	(2) / 99%	0.012	1:10 ⁵	0.1170 μM
FeCl ₃ ·6 H ₂ O	(1) / Analytic	0.012	1:10 ⁵	0.1180 μM
ZnSO ₄ ·7 H ₂ O	(3) / >99.5%	0.080	1:10 ⁸	0.0008 μM
CoCl ₂ ·6 H ₂ O	(1) / Analytic	0.050	1:10 ⁸	0.0005 μM
MnCl ₂ ·4 H ₂ O	(1) / Analytic	0.900	1:10 ⁸	0.0090 μM
Na ₂ MoO ₄ ·2H ₂ O	(2) / ACS	0.030	1:10 ⁸	0.0003 μM
Na ₂ SeO ₃	(2) / ~98%	0.100	1:10 ⁸	0.0010 μM
NiCl ₂ ·6 H ₂ O	(1) / Analytic	0.100	1:10 ⁸	0.0010 μM

Manufacturer Index: (1) Mallinckrodt, (2) Sigma, (3) Fluka

NOTE: To avoid contamination, dry nutrients and metals should be weighed without using spatulas.

1. Prepare nutrient stocks as following using Sigma Ultra grade chemicals:

A. 25mM NaH₂PO₄ (pH 7.5)

- Weigh out 0.15 g NaH₂PO₄·H₂O using dust-free weigh paper
- Transfer into 50 mL Falcon tube filled with about 40 mL Milli-Q water
- Dissolve NaH₂PO₄·H₂O by stirring, inverting tube
- Adjust volume to 50 mL with Milli-Q water
- Filter sterilize with **pre-washed** Millipore Steriflip 0.22 μM Durapore PVDF membrane system (cat: SE179M6); pre-wash with 50 mL milli-Q water and in laminar flow hood, remove from collection tube; reset tube into filtration system
- In a laminar flow hood, apply sterile tube cap.
- Store sterile stock at 4°C

B. 0.4M (NH₄)₂SO₄

- Weigh out 2.64 g (NH₄)₂SO₄ on dust-free weigh paper
- Transfer into 50 mL Falcon tube filled with about 40 mL Milli-Q water
- Dissolve (NH₄)₂SO₄ lume to 50 mL with Milli-Q water
- Filter sterilize with **pre-washed** Millipore Steriflip 0.22 µM Durapore PVDF membrane system (cat: SE179M6); pre-wash with 50 mL milli-Q water and in laminar flow hood, remove from collection tube; reset tube into filtration system
- In a laminar flow hood, apply sterile tube cap.
- Store sterile stock at 4°C

C. 0.5M TAPS (pH 8.0)

- Weigh out 60.75 g TAPS on dust-free weigh paper
- Transfer into 1L glass beaker with stir bar and 300mL Milli-Q water
- Dissolve TAPS by stirring
- Fill remaining volume to 500mL mark with Milli-Q water
- Check pH by pipetting out small amount onto pH paper for rough estimate
- Turn stirrer on and **SLOWLY** add NaOH, check pH by pipetting out small amount onto pH paper for rough estimate periodically
- When pH is around 8.0 use a calibrated pH meter probe that has been Milli-Q water rinsed and soaked several times (add NaOH as needed to obtain pH 7.5)
- Transfer to volumetric flask or graduated cylinder and adjust volume to 500mL
- Filter, using a 500mL sterile filter system (Corning code# 430773) with .22µm nylon filter and polystyrene bottle (pre-wash the filter 2-3 times with 100mL Milli-Q water to remove residual carbon) attached to a vacuum pump (use 1 inch of Hg)
- Store sterile stock at 4°C

D. Working trace metal stock

Use the working trace metal stock from the PRO99 medium recipe, **but for the artificial seawater medium ten times LESS of the stock will be added to the artificial medium compared to PRO99.**

2. Prepare “Turks Island Salt Mix”³ by adding the following:

dissolve each ingredient sequentially in 2/3 final volume of Milli-Q water

Milli-Q water	1.0L
NaCl	28.10g
Mg SO ₄ ·7 H ₂ O	6.90g
MgCl ₂ ·6 H ₂ O	5.49g
KCl	0.67g
CaCl ₂ ·2 H ₂ O	1.47g

- Fill to appropriate volume with Milli-Q water
- Filter-sterilize with a 0.2 µm-cutoff filter (pre-rinsed to remove excess carbon)*

³ The Merck Index (1989) Merck&Co, Inc. (Rahway, N.J.)

3. To Turks Island Salt Mix add the following sterile nutrients stocks:

	<u>1L</u>
NaH ₂ PO ₄ ·H ₂ O	2mL
(NH ₄) ₂ SO ₄	1mL
TAPS	2mL
Working trace metal stock	10µl

4. Store this ASW basal medium at room temperature for up to one month.

5. Prepare 0.6M NaHCO₃ stock fresh on day of medium inoculation (do not store)

- Weigh out 2.52g NaHCO₃ on dust-free weigh paper
- Transfer into 50 mL Falcon tube filled with about 40 mL Milli-Q water
- Dissolve NaHCO₃ by stirring, inverting tube
- Adjust volume to 50 mL with Milli-Q water
- Filter sterilize with **pre-washed** Millipore Steriflip 0.22 µM Durapore PVDF membrane system (cat: SE179M6); pre-wash with 50 mL milli-Q water and in laminar flow hood, remove from collection tube; reset tube into filtration system
- In a laminar flow hood, apply sterile tube cap.
- Store sterile stock at 4°C

6. NaHCO₃ addition to basal medium:

10 mL of 0.6M NaHCO₃ per 1L basal medium

*** AMP-J medium is filter-sterilized; to make "AMP-A", autoclave the Turks Island Salt mix; note that this will increase the concentration of hydrogen peroxide in the medium. Note, AMP-A is not a published recipe (yet)**

III. INSTANT OCEAN MEDIUM

1. Add to 900mL Milli-Q water:

35g Instant Ocean Sea Salts*

2. Bring to 1L with Milli-Q water and stir until dissolution achieved. Autoclave. Allow to cool.

3. Add sequentially to cooled medium:

2.0 mL 0.5M TAPS, pH 8.0

1.6 mL 0.5M NH₄Cl

2.0 mL 0.025M NaH₂PO₄ H₂O

100 µL Trace Metal Working solution

*Use only instant ocean sea salt mix (reef mix has Cu toxic to Pro in it)

Notes- All nutrients prepared as in “procedures for culturing Prochlorococcus” manual. Upon autoclaving white precipitate of excess bicarbonate will come out of solution giving medium a milky white appearance. When aliquotting out IO for culturing purposes agitate solution often to keep precipitate in solution. Unsterilized Instant Ocean has shown evidence of lytic activity (Prochlorococcus-infecting phage) so it is very important to sterilize the solution.

IV. AMP-J AGAR PLATES

Media preparation (1 L AMP ~50 plates)

- A. Noble Agar (0.42% final)
- a. Acid washed polycarbonate bottle w/ minimum volume of 2x MQ-H₂O used (i.e. 500 ml bottle for 250ml agar)
 - i. 4.2g agarose
 - ii. 250 ml MQ-H₂O
 - iii. microwave tindalize:
 1. bring to boil in microwave – **make sure cap is loose!**
 2. agar will violently boil over, watch carefully – stop and restart microwave to get rolling boil.
 3. Carefully remove from microwave, place in hood until agar solidifies or O/N; repeat 3x.
- B. Turks Island Salt Mix Basal Medium
- a. TISM – 750 ml

i. MQ-H ₂ O	750 ml
ii. NaCl	21.075 g
iii. MgSO ₄ •7H ₂ O	5.175 g
iv. MgCl ₂ •6H ₂ O	4.118 g
v. KCl	0.503 g
vi. CaCl ₂ •2H ₂ O	1.103 g
 - b. Autoclave in 1 L acid washed PC bottle
 - c. Prior to last agar tindalization, prewarm Turks Island Salt Mix; place in 37C incubator or remove from autoclave and allow cool to ~37C
- C. Pour AMP plates
- a. Prewarm TISM as noted in step “B. c.”
 - b. Add nutrients to make 1L of 1x AMP (in order shown)

i. Stock NaH ₂ PO ₄	2 ml
ii. Stock (NH ₄) ₂ SO ₄	1 ml (or 1.6 ml NH ₄ Cl stock)
iii. 1 M TAPS pH 8.0	1 ml
iv. TMWS	0.1 ml
v. Fresh 0.6M NaCO ₃ ⁻	10 ml
 - c. gently pour warm agar into 750 ml AMP, swirl gently, pour back and forth from 500 ml bottle to ensure homogeneity of agar.
 - d. Working quickly, pour or use serological pipette to transfer ~20 ml AMP agar to petri dishes. **Note:** Use of serological pipette will make it easy to aspirate any bubbles or imperfections generated in distribution of agar and ensure consistency of agar plates.
- D. Common additions or variations
- a. 0.21% agar – replace 4.2 g Noble agar for 2.1 g. Makes very soft plates.
 - b. Sulfite 20mM final concentration:
 - i. 1M Sulfite made fresh – 1.26 g in 10 ml. filter sterilize
 - ii. 2 ml addition to AMP after bicarbonate addition

V. CULTUREWARE FOR GROWING *PROCHLOROCOCCUS*

Prochlorococcus is highly sensitive to trace metal concentrations, so special care must be given in choosing and preparing media containers and culture vessels. We have had success with the disposable and reusable cultureware, though reusable cultureware requires cleaning with acid and autoclaving (see protocol below).

A. TYPES OF CONTAINERS

Media containers

We use acid washed Nalgene polycarbonate or Teflon bottles.

Culture vessels

We use a range of vessel types. Acid washed polycarbonate or glass bottles or tubes work well, as do disposable polystyrene culture tubes: Falcon 14ml round bottom tube (Becton Dickinson #352051).

B. LIQUID TRANSFER AND MEASUREMENT DEVICES

For small volumes we use sterile disposable plastic pipettes and plastic tips. No acid washing is required. For larger volumes we use acid-washed graduated cylinders (glass or plastic) and volumetric flasks (glass).

C. PROTOCOL FOR ACID WASHING CULTUREWARE

NOTE: under no circumstances should a brush be used to clean the glassware, as this is a significant source of metal contamination. Be sure that there are no rusty objects anywhere in the vicinity of your sink or bench.

- Soak cultureware overnight in a ~2% Micro-90 (Intl. Products Corp.) detergent bath
- Rinse 6X with tap water
- Rinse 6X with deionized water
- Soak overnight in 1N HCl (trace metal free; we use Mallinckrodt #2062)
- Rinse 6X with Milli-Q water (18M Ω water)
- Fill with Milli-Q water and let sit overnight
- Autoclave cultureware filled with Milli-Q water
- Sterilely remove water prior to use

VI. GROWING *PROCHLOROCOCCUS* IN LARGE VOLUME (1-10L)

A. METHOD 1: BUFFERING WITH TAPS AND BICARBONATE:

- Acid wash a 1L or 10L Nalgene polycarbonate container. Add appropriate volume of media
- Supplement medium with TAPS pH 8.0 and sodium bicarbonate
for 1L: 5mM TAPS, 6mM bicarbonate final concentrations
for 10L: 10mM TAPS, 12mM bicarbonate
- Incubate with cap loosened

Note: 0.5 M TAPS and 0.6 M sodium bicarbonate stocks can be made, and should be filter-sterilized with a 0.2 μ m filter. The TAPS can be stored at 4°C for several months, but the bicarbonate should be made fresh the day of inoculation.

B. METHOD 2: PUMPING STERILE AIR OR CO₂ INTO SYSTEM:


- Acid wash a 1L or 10L Nalgene polycarbonate container. Add appropriate volume of media.
- Fix top with filtered air system (see figure, next page)
- Incubate with cap tightened, low rate of airflow (1-4 bubbles/second). A higher airflow rate will cause higher evaporation and the culture may crash as a result

Note: Air source can be an aquarium pump, or compressed air, or CO₂ tank.

System components (for 1L flask):

<u>Part</u>	<u>Vendor</u>	<u>Part number</u>
Nalgene 50 silicone tubing (.187 ID, .312 OD)	VWR	8060-0040
6mm Pyrex glass tubing	Carolina Biol. Suppl.	WW-71-1146
Autoclavable 0.2 μ m filter Acro 37 TF	Pall Gelman	4464
Threaded Teflon cap for 1L flasks (1/4 inch tube diameter)	Savillex	0738-4-2

Bubbling system for *Prochlorococcus*:



VII. PURITY BROTH MEDIA

NOTE: both 1/10X ProAC (rich) and PLAG (minimal) media should be used to routinely test the purity of the axenic *Prochlorococcus* cultures

A. 1/10X ProAC prep

Autoclave separately:

1. 750 ml SSW, filtered with 0.2 μm PC filter
2. 250 ml H₂O + 1.7g AC Difco broth

After autoclaving:

Combine, add nutrients:

- 1.6 ml 0.5M NH₄Cl
- 2.0 ml 0.025 NaHPO₄
- 100 μl trace metal stock
- 0.5 ml 2000X Va vitamins (see below)

B. PLAG prep

Autoclave separately:

1. 750 ml SSW, filtered with 0.2 μm PC filter
2. 250 ml H₂O

After autoclaving:

Combine, add nutrients:

- 1.6 ml 0.5M NH₄Cl
- 2.0 ml 0.025 NaHPO₄
- 100 μl trace metal stock
- 0.5 ml 2000X Va vitamins (see below)
- 2.5 ml Pyruvate (0.05% final)
- 2.5 ml Lactate (0.05% final)
- 2.5 ml Acetate (0.05% final)
- 2.5 ml Glycerol (0.05% final)

Pyruvate, Lactate, Acetate, Glycerol made up as 20% stocks and filter sterilized.

Va vitamin mix (2000X stock):

To prepare stock, dissolve each compound separately, then add together and bring to 1L. Filter sterilize, store in the dark at 4°C.

Vitamin	Final Stock Conc. (2000X)
Thiamin-HCl	0.2g/L
Biotin	0.001 g/L
Vitamin B ₁₂	0.001 g/L
Folic acid	0.002 g/L
PABA (para-aminobenzoic acid)	0.01 g/L
Nicotinic acid	0.1 g/L
Inositol	1.0 g/L
Calcium pantothenate	0.2 g/L
Pyridoxine-HCl	0.1 g/L

Reference: Waterbury and Willey. Methods in Enzymology 167: 100-105.

C. Purity assay

Routinely add 1 ml of axenic *Pro* cultures to 5-8ml 1/10 ProAC and PLAG medium.

Incubate at room temperature in dark and monitor for growth over 1+ month
(note: *Pro* will grow in PLAG in the light, so it is important to keep these tubes in the dark)

VIII. CRYOPRESERVATION PROTOCOL

A. REUSABLE FROZEN STOCK (CONCENTRATED CELLS):

1. Preparation of (50X) concentrated cells

- Grow cultures in light level approx. $16 \mu\text{Ein}/\text{m}^2/\text{sec}$, continuous light, and $20\text{-}22^\circ\text{C}$ temperature [if large volume (500ml+) add 6 mM sodium bicarbonate and 5mM HEPES pH 7.5].
- At mid log to early stationary phase, transfer volume of cells (50mL) into sterile centrifuge tubes/bottles.
- Spin at 8000RPM for 15 minutes at $20\text{-}22^\circ\text{C}$.
- Carefully decant liquid into waste container.
- Resuspend pellet in 1ml fresh medium (same medium used for growth).
- Withdraw concentrated cells (1mL) and transfer into 1.2mL cryogenic vial (Nalgene #5000-0012).
- Add $75 \mu\text{l}$ DMSO (7.5% final concentration) to cryogenic vial and invert several times.
- Immediately place into liquid nitrogen or into -150°C freezer.
- Note: recovery from liquid nitrogen is faster, but storage in liquid runs the risk of cross-contamination from other vials.

2. Recovery of cells in fresh medium

- Reduce ambient light in work area as much as possible.
- Remove cryogenic vials one at a time from liquid nitrogen or freezer and transport vials under cover to laminar flow hood as quickly as possible.
- With a sterile toothpick, immediately scrape up a small amount of the frozen cells and place into 25 mL sterile medium.
- Note: autoclaved wooden toothpicks work fine, as do ethanol-sterilized plastic toothpicks.
- Once scrape is completed return cryogenic vial to liquid nitrogen/ -150°C freezer.
- Incubate cultures at $20\text{-}22^\circ\text{C}$ under low continuous light ($\sim 16 \mu\text{Ein}/\text{m}^2/\text{sec}$).
- Note: **Recovery under low light is absolutely critical.**
- Monitor growth with daily fluorescence measurements. Growth can be detected within one week after thawing, but may take as long as 3-4 weeks.

B. SINGLE USE FROZEN STOCK (UNCONCENTRATED CELLS): SUITABLE FOR SHIPMENTS ON DRY ICE

1. Preparation of unconcentrated cells

- Grow cultures in light level approx. $16 \mu\text{Ein}/\text{m}^2/\text{sec}$, continuous light, and $20\text{-}22^\circ\text{C}$ temperature (if large volume add 6 mM sodium bicarbonate and 5mM TAPS pH 8.0).
- At mid log to early stationary phase, transfer 1 ml cells into a 1.2mL cryogenic vial (Nalgene #5000-0012).
- Add 75 μl DMSO (7.5% final concentration) to cryogenic vial and invert several times.
- Immediately place into liquid nitrogen or into -150°C freezer.
- Note: Recovery from liquid nitrogen is faster, but storage in liquid runs the risk of cross-contamination from other vials.

2. Recovery of cells in fresh medium

- Reduce ambient light in work area as much as possible.
- Remove cryogenic vials from liquid nitrogen or freezer and transfer to 37°C bath.
- Partially thaw the cells: once half of liquid has thawed and small ice chunk remains
- Remove vials and transport to laminar flow hood.
- Transfer $2 \times 500 \mu\text{l}$ cells into 2 culture tubes containing 25 mL sterile medium.
- Incubate cultures at $20\text{-}22^\circ\text{C}$ under low continuous light ($\sim 16 \mu\text{Ein}/\text{m}^2/\text{sec}$).
- Note: **Recovery under low light is absolutely critical.**
- Monitor growth with daily fluorescence measurements. Growth can be detected within one week after thawing, but may take as long as 3-4 weeks.

Note: Both recovery methods will work for concentrated cells, and may work for unconcentrated cells. Significantly faster recovery of unconcentrated cells occurs when the entire contents of the vial are used for inoculation.

Note: Both PRO99 and AMP-J media are suitable for preparing cells for freezing and for recovery after thawing.